



# **EUV source metrology using transmissive and diffractive optics**

Muharrem Bayraktar  
[m.bayraktar@utwente.nl](mailto:m.bayraktar@utwente.nl)



EUV Source Workshop  
23 October 2024

# Motivation




- Light outside the 13.5 nm ( $\pm 1\%$ ) band, spans a broad spectrum (DUV, VIS, NIR)
- Can be reflected from the optics and cause contrast loss, and can be absorbed by H<sub>2</sub>
- Need for broadband (both in-band and out-of-band) spectral and spatial monitoring

# In-band metrology: Flying Circus approach






- Absolutely calibrated measurement in the pass-band of the mirror, <5% calibration uncertainty



- Enabled a through comparison of early EUV sources in mid-2000s,
- Varieties are in common use nowadays

F. Bijkerk, et.al. EUV Sources for Lithography, 721-734, (2006)

# In-band metrology: Anomalous transmission filters





- Based on anomalous transmission (Borrmann) effect
- Can be coated on photodiodes
- 74 bilayers of Ni/Si – transmission 1%, bandwidth FWHM 0.23nm (1.7%) at 13.5nm wavelength
- Passband can be tuned by changing the incidence angle
- Intermixing plays a significant role in the final response, can be improved

J. L. P. Barreaux, et.al. Opt. Exp. 25(3), 1993 (2017)

# Broadband spectroscopy: Transmission grating spectrometer

- ✓ Grating chip developed in MESA+ Nanolab
- ✓ 3x7 grating matrix




Low resolution gratings (500 lines/mm) for VIS range  
High res. gratings (up to 10.000 lines/mm) for EUV range






# XUV transmission grating spectrometer

- ✓ Broadband coverage
- ✓ Compact design ( $\sim 20 \times 25 \times 30 \text{ cm}^3$ )




- ✓ Computer controlled positioning
- ✓ Straightforward alignment



M. Bayraktar, et.al. NEVAC Blad, 54, 14-19 (2016).



# Scaling of emission in an industrial EUV source



Specific charge states can be resolved clearly in the short wavelength-side of the 13.5nm peak

F. Torretti, et.al. J. Phys. D: Appl. Phys. 53, 055204 (2020)

# Absolutely calibrated spectra in the 5.5-265 nm range



- ✓ Strong higher diffractions orders of the 13.5nm peak are suppressed by the set of filters
- ✓ SiC 4-bounce filter enables to measure the otherwise inaccessible 40-115nm range


Z. Bouza, et.al. AIP Adv. 11, 125003 (2021)



# Broadband spectroscopy and imaging using zone-plates

*Y. Mostafa, Z. Bouza, et. al. Opt. Lett. 48, 4316-4319 (2023)*

- Pinhole imaging with foil filters – broad spectral coverage but limited spatial resolution
- Imaging using multilayer mirrors – decent spatial resolution but limited spectral coverage
- **Combination of transmission gratings with dispersion matched zone-plates – broadband spectral coverage + decent spatial resolution**




- Spectroscopy on x-axis
- Imaging on y-axis
- Light focusing using zone-plates is strongly wavelength dependent
- How to have good focusing for a broad wavelength range (e.g. 5-80 nm)?

# Dispersion-matched zone-plates enable broadband spectroscopy and imaging


1D zone-plate designed to focus a single wavelength

$$r_{i,n} = \sqrt{n\lambda_i \left( f + \frac{n\lambda_i}{4} \right)}$$

1D tapered zone-plate designed to focus 5-80 nm wavelength range




- Radius  $r_{i,n}$  of each zone  $n$  is dependent on the wavelength  $\lambda_i$
- As the wavelength increases, the zone radius increases
- Resulting in tapered zone plates



Y. Mostafa, et. al. Opt. Lett. 48, 4316 (2023)


# Experimental implementation



- Nd-YAG laser at 1  $\mu\text{m}$
- Sn solid target
- Beam size: 60  $\mu\text{m}$  FWHM
- Magnification: 1.9
- Design spatial resolution: 10  $\mu\text{m}$
- Spectral resolution: 0.8 nm @13.5nm
- Better spectral and spatial resolutions possible by setup configuration
  - Spectral resolution: < 0.1 nm
  - Spatial resolution: < 1  $\mu\text{m}$

Y. Mostafa, et. al. Opt. Lett. 48, 4316 (2023)


# Results – experiment vs simulation



- Spectrum without zone-plate, as a reference
- Spectral resolution 0.8 nm
- Prominent 13.5 nm peak is visible
- Experimental recording with zone-plate inserted
- Horizontal band at the center corresponds to the focused plasma light
- Zero- and other zone-plate orders are visible at both sides of the plasma light
- Wave propagation simulation, assuming zone-plate is illuminated with a plane-wave, yielding the theoretical design resolution (10  $\mu\text{m}$ ) of the zone-plate


Y. Mostafa, et. al. Opt. Lett. 48, 4316 (2023)

# Broadband imaging and spectroscopy in the 5-80 nm range using zone-plates





# Broadband imaging and spectroscopy in the 5-80 nm range using zone-plates



- Broadband imaging and spectroscopy in the 5-80 nm range is possible
- Normalized image shows a widening central band, indicating increasing plasma size with increasing wavelength as anticipated

# Conclusions

- In-band metrology using Flying Circus and Anomalous transmission filters approaches
- Broadband metrology using transmission gratings
  - Broadband spectroscopy in the EUV-NIR range, compact setup with good spectral resolution
  - Measurements using HHG, prototype EUV light source (ASML) and research type LPP sources (ARCNL)
  - Enables to observe distinct emission lines of the tin plasma
- Dispersion-matched zone-plates enable broadband spectroscopy and imaging
  - 0.8 nm spectral resolution and 10  $\mu\text{m}$  design spatial resolution
  - Broad spectral coverage 5-80 nm
  - Better spectral ( $<0.1$  nm) and spatial ( $<1$   $\mu\text{m}$ ) resolutions possible
  - Wave propagation simulations verify the spatial resolution and the proof-of-principle of the tool
  - Laser-produced-plasma investigated with the tool show increasing plasma size with increase wavelength, as expected

# Acknowledgements

## Spectrometer Team

Dr. James Byers (UTwente, NL)

Dr. Bert Bastiaens (UTwente, NL)

Dr. Boris Vratzov (NT&D, DE)

Caspar Bruineman (Scientec, NL)

Prof. Dr. Fred Bijkerk (UTwente, NL)

Prof. Dr. Marcelo Ackermann (UTwente, NL)

## Collaboration with ASML

Dr. Fei Liu

Dr. Mark van de Kerkhof

ASML Veldhoven Source Performance Team

Dr. Igor Fomenkov and Cymer Team

## Collaboration with ARCNL, NL

Dr. Zoi Bouza

Dr. Ievgeniia Babenko

Dr. Yahia Mostafa

Dr. Francesco Torretti

Dr. Ruben Schupp


Dr. Lars Behnke

Prof. Dr. Wim Ubachs

Prof. Dr. Oscar Versolato

# Appendix

# Resolution test in a HHG source




- Design resolution 0.09 nm at 21 nm (37<sup>th</sup> harmonic)
- Measured resolution  $< 0.13$  nm<sup>[5]</sup>
- Can be configured for 0.05 nm resolution at 13.5 nm

[5] J. Goh et.al. Opt. Exp. 23, 4421 (2015).




# Demonstration in the DUV, VIS, NIR range



- Low but detectable DUV-NIR intensity
- Distinct peaks and broad features in the spectra

# Assigning emission lines



- Peaks originate from higher diffraction orders, line and continuum emission
- Higher diffraction orders can be filtered out by spectral windows (eg. LiF, MgF<sub>2</sub>, UVFS...)
- Sn<sup>+2</sup>, Sn<sup>+3</sup> and H lines correlate with observed spectra